

**Product Name :**  
Ambu Type Bag (Artificial Resuscitator), (Autoclavable)  
Superior, Adult, Black Rubber

**Product Code :**  
NMANI10039



**Description :**

Ambu Type Bag (Artificial Resuscitator), (Autoclavable) Superior, Adult, Black Rubber Black rubber. Optional : Reservoir bag - AN403 (Order separately). CE certified

**Technical Specification :**

Rubber Ambu Bag (Resuscitator) kit consists of manual Resuscitator, Tubing, Reservoir Bag & Mask.

It provides excellent bag re-expansion property. Pressure limiting valve ensures crucial adjustments for enhanced performance.

It is available in two sizes : Adult & Child.