

Product Name :
Ambu Type Bag, (Artificial Resuscitator), Child Black
Rubber with special valve

Product Code :
NMANI10040



Description :

Ambu Type Bag, (Artificial Resuscitator), Child Black Rubber with special valve Black rubber. With face mask size. 00 (infant) & 1 (round). CE certified

Technical Specification :

Rubber Ambu Bag (Resuscitator) kit consists of manual Resuscitator, Tubing, Reservoir Bag & Mask.

It provides excellent bag re-expansion property. Pressure limiting valve ensures crucial adjustments for enhanced performance.

It is available in two sizes : Adult & Child.