

<b>Product Name :</b> Ambu Type Bag, (Artificial Resuscitator), Child Black Rubber with special valve - Autoclavable	<b>Product Code :</b> NMANI10041
---	-------------------------------------



<b>Description :</b>
Ambu Type Bag, (Artificial Resuscitator), Child Black Rubber with special valve - Autoclavable Black rubber. With face mask size. 2. Round shape. CE certified
<b>Technical Specification :</b>
Rubber Ambu Bag (Resuscitator) kit consists of manual Resuscitator, Tubing, Reservoir Bag & Mask.  It provides excellent bag re-expansion property. Pressure limiting valve ensures crucial adjustments for enhanced performance.  It is available in two sizes : Adult & Child.