

Product Name : Ambu Type Bag, (Artificial Resuscitator), Child Black Rubber with special valve - Autoclavable	Product Code : NMANI10041
---	-------------------------------------



Description :
Ambu Type Bag, (Artificial Resuscitator), Child Black Rubber with special valve - Autoclavable Black rubber. With face mask size. 2. Round shape. CE certified
Technical Specification :
Rubber Ambu Bag (Resuscitator) kit consists of manual Resuscitator, Tubing, Reservoir Bag & Mask. It provides excellent bag re-expansion property. Pressure limiting valve ensures crucial adjustments for enhanced performance. It is available in two sizes : Adult & Child.