

Product Name :
Oxygen Cylinder, Steel, Empty

Product Code :
NMANI10052



Description :

Oxygen Cylinder, Steel, Empty

Technical Specification :

Gas cylinders are made of Seamless Steel tubes, in Manganese steel and Chrome Molybdenum Steel.

The technology used for manufacturing medical gas cylinders is hot spinning of seamless steel tubes for closing both ends and subsequently followed by other finishing processes.

The Equipments for manufacturing medical gas cylinders have been sourced from suppliers who are world leaders in their respective fields.