

Product Name :
Anesthesia machine with LED Screen**Product Code :**
NMANI10054**Description :**

Anesthesia machine with LED Screen With LED screen. Mounted on 4 castors. 2 front castors with brakes. With 2 big drawers . Handle for easy maneuvering Epoxy painted high quality mild steel frame. Equipped with Anesthesia Vaporizer. Agent type: Isoflurane, Enflurane, Halothane. Classical CO₂ circle absorber with 1 canister With necessary monitoring and alarm function, to ensure the operational safety

Technical Specification :

Suitable for adults and pediatric patients
Equipped with Large LCD screen, various monitoring waveform and parameters indifferent colors
Pneumatically driven electronically controlled
Integrated CO₂ absorption system equipped with NO WATER VAPOR COAGULATION sensor and fast release APL valve
6-tube of flow meters for 3-gas sources meet long operation requirement
Spirometry loops
Equipped with Anesthesia Vaporizer.

