

**Product Name :**  
Anesthesia machine integrated for Pediatric and Adult**Product Code :**  
NMANI10058**Description :**

Anesthesia machine integrated for Pediatric and Adult With 10.4 TFT color screen Mounted on 4 castors. 2 front castors with brakes. With 1 safety lock and 3 big drawers. Handle for easy maneuvering. Aluminium alloy rail frame. With durable and non-corrosive stainless steel platform. Equipped with high precise Anesthesia Vaporizer (with temperature, flow and pressure compensation function). The agent output is continuous, stable & precise The agent type: Isoflurane, Enflurane, Halothane, Sevoflurane. With Volume & pressure ventilation modes Classical CO2 circle absorber with 2 canisters for low resistance and minimum gas leakage Fresh gas outlet is in front of frame for the user to easily connect Jackson-Rees circuit. With upgraded with functions, eg, CO2 Monitor, Multi-parameter monitor, Anesthesia depth monitor, etc Working interfaces for adult and pediatric are separated; the doctor can choose interface according to patient type Optional: O2, N22 Yoke device. When main gas supply is fails, The backup gas supply starts.

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website: SITE,