

Product Name :
Rack Steam Sterilizer, Single Rack**Product Code :**
NMALV10016**Description :**

Rack Steam Sterilizer, Single Rack Steam sterilizer, pressure type, single rack, used to sterilize reusable: syringes and needles. Aluminium vessel with high pressure seal suitable for carrying out sterilisation using superheated steam. Steam sterilizer designed to operate on external heat source (non-electric) at +121½C over a 20 minute period. Components: Aluminium vessel with heat resistant handles. Aluminium cover with heat resistant handles , rubber sealing gasket, pressure valve & safety plug. 1 aluminium rack with 1 lid, for 36 syringes (4 syringes 5ml, 16 syringes 2 ml, 16 syringes 1 ml) and 40 needles. Sterilization capacity: 5.5 litres
Supplied with: 1 x sterilization indicator sticker (pack of 100), 2 x spare safety plugs, 3 x spare sealing gaskets, 1 x plastic bowl for washing syringes, 1 x carrying bag, 1 x user manual in English.

Technical Specification :

Made of Aluminium
For placing of syringes inside autoclave
Holds 36 syringes (4 syringes 5 ml, 16 syringes 2 ml, 16 syringes 1 ml) and 40 needles
With Lid

