

**Product Name :**  
Accessories for Rack Sterilizer**Product Code :**  
NMALV10018**Description :**

Accessories for Rack Sterilizer

**Technical Specification :**

Three different types of stainless steel containers of different heights (large, medium and small) are available for handling different loads.

Stacking shorter containers can increase load capacity (two medium-sized and three small-sized containers can fit in a given chamber).

A mesh wire basket is suitable for most loads, including glassware and media.

A pail or a mesh wire basket with a solid bottom is used for agar and to contain potential spillage from a torn biohazard bag.

A drum has a hinged lid with a handle. Depending on the load, it serves as a closed (but not watertight) trash container or as a carrying and storage case.

The side and bottom of the drum are perforated for better steam penetration during sterilization.

Sliding panels cover these perforations prior to and after sterilization.

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,