

Product Name :
Autoclaves Pressure Steam Sterilizers (Aluminum) -
Economy Mirror / Polished Finish

Product Code :
NMALV10019



Description :

Autoclaves Pressure Steam Sterilizers (Aluminum) - Economy Mirror / Polished Finish Construction: These autoclaves are made of 10 gauge aluminum. The drum is deep drawn & the ring on top is riveted. The lid is made by casting. Locking of lid by 6 nuts with rubber gasket seal in between (not metal to metal seal). Two safety valves, pressure gauge, flexible metal rubbing to exhaust air. Includes outer stand, inner aluminum container, inner tripod stand.

Technical Specification :

Type: Pressure Cooker Type/Laboratory bench autoclave(N-Type). Dimensions: 12" Dia X 12" Height.
Capacity: 23 ltrs.

Equipped with pressure gauge, steam safety valve, pressure regulation device, steam release valve, vacuum release valve with heat resistance lid handle and body handle, outer stand.

Ideal Use: To sterilize medical, equipment, glassware, surgical, equipment, medical wastes, culture media, autoclave containers, plastic tubes, and pipette tips.

