

**Product Name :**  
Autoclave 11 Litre Electric

**Product Code :**  
NMALV10024



**Description :**

Autoclave 11 Litre Electric Stainless Steel 304 seamless construction. Dimension : 220mm Dia x 300mm.Ht. (with Lid) Note: Height and capacity adjustable as per customers requirement

**Technical Specification :**

All materials & components of food grade quality.

Precisely tested pressure regulator and precision pressure gauge.

Dressing Drum Stainless Steel SEAMLESS

Gasket Release System (GRS) and safety valve that automatically defuses pressure to safe level.

Complete with Double safety valve, Steam valve,& water Release Cock Pressure Gauge