

**Product Name :**  
Light Weight Disinfection Tray System**Product Code :**  
NMALV10030**Description :**

Light Weight Disinfection Tray System This tray system is a nice sterilizing / disinfecting / soaking tray that is essential for sterilization of certain instruments & equipment which cannot be sterilized under steam & pressure in Autoclaves. These can be also used as a back up to your other sterilizing system. The tray includes Outer Tray with lid & a Lift-out drain basket with handles to facilitate lifting of instruments in and out of the soaking solution. Lid protects against spillage and minimizes vaporization Dimension: Outer tray measures 845mm L x 225mmWx 150mm H Volume Capacity: 10 liters appx, [NON AUTOCLAVABLE]

**Technical Specification :**

This tray system is the nice sterilizing / disinfecting / soaking tray that is essential for sterilization of certain instruments & equipment which cannot be sterilized under steam & pressure in Autoclaves.

These can be also used as a back up to your other sterilizing system.

The tray includes Outer Tray with lid & a Lift-out drain baskets with handles to facilitate lifting of instruments in and out of the soaking solution.

Lid protect against spillage and minimize vapors.

---

## Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,

**NAUGRA<sup>®</sup>**