

**Product Name :**  
Double Wall Semi Automatic Autoclave**Product Code :**  
NMALV10032**Description :**

Double Wall Semi Automatic Autoclave

**Technical Specification :**

Inner & Outer Chamber, Constructed from S.S 304 Grade  
Bucket / Basket Made of Stainless Steel  
Water Level indicator  
Lid Closing by Radial Locking System  
Sealing of lid by Neoprene Rubber Basket  
Lid is fitted with Pressure Gauge, Spring Loaded Safety Valve  
Operating range from 15psi to 20psi pressure  
Sterilizing temperature upto 134 c  
Voltage:220 VoltsAC (50 Hz)  
Tested upto 2.5 times of required pressure  
Fully automatic models requiring no human invention till a full sterilization cycle is over.  
Fitted with a (Programmable Logic Controller) LCD Display with respect to time and temperature  
Enabling user automatic purging of stale air.  
sterilizing holds time system for sterilizing period and cycle end.  
Automatic exhaust of Pressure  
Audio alarm/ buzzer is raised at the end of the cycle on.  
All internal joints are argon arc welded, ground and polished to give crevice free internals.  
Gasket is made from imported silicon rubber.  
Drain valve facilitates easy cleaning of the chamber.

---

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

**Website:**Medical Hospital Equipments Manufacturers,