

Product Name :
Front Loading Autoclave**Product Code :**
NMALV10033**Description :**

Front Loading Autoclave Class-N Front Loading Autoclave. With thermal vacuum, fully automatic with dry cycle. 15 Lt. capacity with Micro Processor control. Digital Display of Pressure & temperature Double locker on door safety feature A pressurized steam autoclave, with thermal air ventilation & drying cycle, handles N type instruments. Automatic diagnostic & alarm system is installed

Technical Specification :

A front-loading autoclave are the best option available for research, education, and quality assurance facilities.
Table top front loading autoclave.
Fully automatic micro-processor based control.
Sterilization chamber made up of thick stainless steel sheet, tested for high pressure.
Digital timer both for wet & dry cycle.
Auto drain of water to reservoir tank & also condensation.
Alarm after the total command is completed.
Cycle time 10-15 mins.

