

Product Name :
Hospital Blood Bank Refrigerator**Product Code :**
NMBBE10001**Description :**

Hospital Blood Bank Refrigerator Tender Suppliers Africa and Hospital Blood Bank Refrigerator Manufacturers Africa Kenya, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Tunisia, South Sudan, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, Benin, Burundi. Medical Laboratory Instruments, Hospital Lab Instruments Equipments, Hospital Laboratory Instruments Tools, Medical Lab Instruments Tool Kits, Hospital Blood Bank Refrigerator Manufacturers, and Hospital Laboratory Instruments Suppliers from India, Hospital Blood Bank Refrigerator China, Hospital Blood Bank Refrigerator Kenya.

Hospital Blood Bank Refrigerator Tender Suppliers, Guadeloupe (Basse-Terre), Guatemala (Guatemala City), Guyana (Georgetown), Haiti (Port-au-Prince), Honduras (Tegucigalpa), Jamaica (Kingston), Martinique (Fort-de-France), Mexico (Mexico City), Montserrat (Brades), Nicaragua (Managua), Panama (Panama City), Paraguay (Asunción), Peru (Lima), Puerto Rico (San Juan), Saint Barthelemy (Gustavia), Trinidad and Tobago (Port-of-Spain), Turks and Caicos Islands (Cockburn Town), USA (Washington DC, America), Uruguay (Montevideo), Venezuela (Caracas), Virgin Islands (Charlotte Amalie), Anguilla (The Valley), Antigua and Barbuda (Saint Johns), Argentina (Buenos Aires), Aruba (Oranjestad), Bahamas (Nassau), Barbados (Bridgetown), Belize (Belmopan), Bermuda (Hamilton), Bolivia (Sucre), Brazil (Brasilia), British Virgin Islands (Road Town), Canada (Ottawa), Cayman Islands (George Town), Chile (Santiago), Colombia (Bogota), Costa Rica (San Jose), Cuba (Havana), Curacao (Willemstad), Dominica (Roseau), Hospital Blood Bank Refrigerator Tender Suppliers, Dominican Republic (Santo Domingo), Ecuador (Quito), El Salvador (San Salvador), Falkland Islands (Stanley), French Guiana (Cayenne), Greenland (Nuuk), Grenada (Saint Georges), Saint Kitts and Nevis (Basseterre), Saint Lucia (Castries), Saint Martin (Marigot), Saint Pierre and Miquelon (Saint-Pierre), Saint Vincent and the

Grenadines (Kingstown), Sint Maarten (Philipsburg), South Georgia and South Sandwich Islands (King Edward Point), Suriname (Paramaribo)

Technical Specification :

Hospital Blood Bank Refrigerator - Blood Bank Refrigerator Microprocessor based controlling with temperature reading is inbuilt in blood bank refrigerators. Blood bank refrigerators have PT-100 sensor for accurate temperature sensing. High and low limit alarming in blood bank refrigerators provides bag safety. Chart recorder with rechargeable battery backup in blood bank refrigerators records temperature during power failure also. Options of ink less/ strip chart recorder are available with blood bank refrigerator. Electric supply for blood bank refrigerator is 220V, 1N, 50HZ. Blood bank refrigerator has inbuilt voltage stabilizer.

Our blood bank refrigerator are available in different capacities; 50 blood bags,110 blood bags,200 blood bags,250 blood bags,360 blood bags. Hospital Blood Bank Refrigerator, Blood Bank Equipments Exporters, Blood Bank Equipments Manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,