

Product Name :
Portable Autoclave - Semi Automatic**Product Code :**
NMALV10039**Description :**

Portable Autoclave - Semi Automatic SS Seamless die pressed deep drawn chamber μ p based PID controller with dual display for Temp. & Time Fully SS construction Dual safety valves Microprocessor Digital temperature controller cum timer Over heat auto shut down at 135 \pm 1/2 $^{\circ}$ C Dual Display Unique shape heater of outer SS tubing Food grade silicon rubber gasket for better durability Specially designed heater for optimum chamber usage Calibrated temperature-cum-pressure gauge Colour coded monitoring gauge SS Safety valve Vacuum breaker to remove air pockets in the chamber

Technical Specification :

Providing a very compact and simple autoclave for the usage of sterilization.
Top cover clamped to the lower tank by wing nut closing system
The unit has a spring-loaded safety valve for pressure adjustment and a dead weight safety valve Supplied with electric fitting, cord, and plug.
Top cover seals the autoclave with neoprene rubber gasket