

**Product Name :**  
Portable Autoclave - Semi Automatic**Product Code :**  
NMALV10039**Description :**

Portable Autoclave - Semi Automatic SS Seamless die pressed deep drawn chamber  $\mu$ p based PID controller with dual display for Temp. & Time Fully SS construction Dual safety valves Microprocessor Digital temperature controller cum timer Over heat auto shut down at  $135\pm\frac{1}{2}^{\circ}\text{C}$  Dual Display Unique shape heater of outer SS tubing Food grade silicon rubber gasket for better durability Specially designed heater for optimum chamber usage Calibrated temperature-cum-pressure gauge Colour coded monitoring gauge SS Safety valve Vacuum breaker to remove air pockets in the chamber

**Technical Specification :**

Providing a very compact and simple autoclave for the usage of sterilization.  
Top cover clamped to the lower tank by wing nut closing system  
The unit has a spring-loaded safety valve for pressure adjustment and a dead weight safety valve Supplied with electric fitting, cord, and plug.  
Top cover seals the autoclave with neoprene rubber gasket