

Product Name :
Horizontal Sterilizer**Product Code :**
NMALV10040**Description :**

Horizontal Sterilizer GMP Compliant Semi automatic operation Pressure switch control High pressure high vacuum option Low water level cutoff Online validation option

Technical Specification :

Triple walled construction inner chamber is made of high grade stainless steel sheet of SS-304 grade.

The outer body is also made of thick stainless steel sheet. The boiler is made of S. Steel.

The lid is made of thick steel sheet plate.

The unit is duly fitted with highly efficient M.S. Painted Ring & radial locking system duly chrome plated.

A three way valve or a multi port Valve made of stainless steel is installed for controlling the sterilization cycle.

The pressure is adjustable from 5 psi to 20 psi, pressure is controlled by automatic pressure control switch .

Fitted with water level indicator and vacuum breaker.

The unit is fitted with temperature gauge, vacuum breaker and steam trap.

The whole unit is mounted on a robust tubular stand made of MS duly painted with attractive stove enamel.

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,