

Product Name :
Fully Automatic Autoclaves**Product Code :**
NMALV10042**Description :**

Fully Automatic Autoclaves This is the fully automatic microchip controlled autoclave. It's a unique vessel that generates steam under pressure and temperature between 120½ to 125½ Celsius. On reaching the right temperature the instruments are subject to a saturated steam environment for a pre specified time which results in safe and guaranteed sterilization. It is sturdily constructed in rust-proof, heavy gauge aluminium sheet. "This autoclave" offers a single step operation for quick and easy use and light sequences shows the cycle's progress. This autoclave offers the most cost effective solution for Steam sterilization Light, Compact, Portable, easy to use top loading autoclave. All Aluminum deep drawn seamless construction. This autoclave offers single step control for quick and easy use. Fast fully automatic microchip controlled sterilization cycle at 121½ C Equipped with double "Over Pressure" Protection system Extended life lid gasket and heating element Color-coded Pressure gauge show internal chamber condition during the cycle

Technical Specification :

In this autoclave we provide automatic purging system through solenoid valve with automatic steam release after sterilization

Supplied with PID Controller with timer to set time for sterilization

Automatic release of pressure after completion of cycle with set timer

Low water level cut of system with alarm for safety of heaters

Specially designed heater for low power consumption

Equipped with 2 Safety Valve, Pressure gauge on top with lid

Max working pressure set on Valve 15 PSI with setting of controller at 121C

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,