

**Product Name :**  
Laboratory Round Heating Plate**Product Code :**  
NMALV10050**Description :**

Laboratory Round Heating Plate The Body is made of P.C.R.C. sheet painted with attractive stove enamel. The top of plate is 20.5 cm dia. fitted with an indicating lamp. This type of heating plate is very useful for fast heating upto 500 $\frac{1}{2}$ C. To work on 220/230 Volts A.C

**Technical Specification :**

These hot plates are manufactured using high grade raw material and consists of cast iron table top, frame made of P.C.R.C. sheet coated powder paint and energy controller to regulate the temperature from ambient up to 350 degree C.

Moreover, the table top of the unit consists of heating elements of 1.5 kW rating and is round in shape with diameter of 20 cm.