

Product Name :
Micro Mist Aerosol Generator**Product Code :**
NMALV10014**Description :**

Aerosol generators to generate droplet size in sub-microns that ensures 100% efficacy. An autoclavable tank covers up to 5000 cubic feet volume in a single fill. The generators are fitted with UL & CE approved imported motors & timer-switch that cuts-off precisely at the pre-set volume of liquid charged. This facilitates unattended operation.

Technical Specification :

Best grade Stainless Steel tank and hardware
All materials compatible with disinfectants in the acidic as well as alkaline range (pH 1.5 to 14)
Anti-skid gasket for stable operation and avoid physical damage S
State of the art Stainless Steel Nozzle for maximum coverage and efficacy
Medical grade silicon tubing

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [SITE](#)