


Product Name :
Hospital Blood Collection Monitor**Product Code :**
NMBBE10004**Description :**

Hospital Blood Collection Monitors - Blood Collection Monitor is a compact instrument to provide smooth and gentle rocking for homogeneous mixing with anticoagulant without clot formation of blood cells during collection of blood from a donor.

- Volume can be set in 1 ml. increments.
 - Provision of pausing collection and of changing programmed volume during pause.
 - Micro-controller based program.
 - Display of weight and volume.
 - Auto tare facility to accurate for the weight of the bag.
 - Motor activated clamping at the end of the collection.
 - Audio visual alarm to alert in case of any abnormal condition.
 - Auto calibration.
 - Highly user-friendly system with LCD messages.
 - Display of the time taken for collection.
 - Diagnostic mode for technical analysis.
- Hospital Blood Collection Monitor, Blood Bank Equipments Exporters,
Blood Bank Equipments Manufacturers

Website: [SITE](#)

NAUGRA[®]

Medical Hospital Equipments Manufacturers