


**Product Name :**  
Hospital Blood Collection Monitor**Product Code :**  
NMBBE10004**Description :**

Hospital Blood Collection Monitors - Blood Collection Monitor is a compact instrument to provide smooth and gentle rocking for homogeneous mixing with anticoagulant without clot formation of blood cells during collection of blood from a donor.

- Volume can be set in 1 ml. increments.
  - Provision of pausing collection and of changing programmed volume during pause.
  - Micro-controller based program.
  - Display of weight and volume.
  - Auto tare facility to accurate for the weight of the bag.
  - Motor activated clamping at the end of the collection.
  - Audio visual alarm to alert in case of any abnormal condition.
  - Auto calibration.
  - Highly user-friendly system with LCD messages.
  - Display of the time taken for collection.
  - Diagnostic mode for technical analysis.
- Hospital Blood Collection Monitor, Blood Bank Equipments Exporters,  
Blood Bank Equipments Manufacturers

---

Website: Medical Hospital Equipments Manufacturers,

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**