

Product Name :
Weighing Scale, Standing, Adult, Digital with Height Rod

Product Code :
HW39



Description :

Weighing Scale, Standing, Adult, Digital with Height Rod Capacity: 300 Kg X 0.1Kg Accuracy: ±150g Weight Unit: Kg/lb LCD Display:1.0" LCD digital Display Dimension: Platform Size: 310 x 310 x 54mm, Column Height:900mm Key Function: Hold/BMI, Unit/Weight, On/Zero/Off, Tare Power Supply: 1.5V 300mA adapter 1.5V AAA battery *6 Operation Temperature: 5°C / 35°C Height Rod: 60-205cm / 6ft 9inch