

**Product Name :**  
Bathroom Scale electronics flat with tare parent model**Product Code :**  
HW45**Description :**

Bathroom Scale electronics flat with tare parent model Maximum Capacity: 300 kg / 660 lbs. Graduation: 0.1 kg / 0.2 lb. Accuracy:  $\pm 150g$  Resolution: 1:3000 Tare Range: Full Range Repeatability(STD dev):  $\pm 1d$  Linearity:  $\pm 3d$  Unit of Measure: Kg/lb Key Function: Hold/Unit/Weight, On/Zero/Off, Tare Stabilisation Time:1-2Sec Operating Temperature: 5°C-35°C or 13° F -95°F Power Supply: 9V battery or AC adaptor Calibration: Automatic External Calibration Mass: 300KG Indicator Display:1.0"LCD digital display Dimension(w x d x h): 15.5 x12.5x3.3 inch 39.5 x 32 x 8.5 cm Gross Weight : 7.7lb/3.5 Kg

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**Website: [SITE](#)