

Product Name :
Vacuum Extractor - Manual

Product Code :
Gynecological0010



Description :

Vacuum Extractor - Manual Vacuum Extractor Set, manual operated. Complete with 3 stainless steel cups (40, 50, 60mm) & 2 Silicon Cups (without valves) Autoclavable Silicone pressure tube. Traction handle. Suction pump Glass jar fitted between M.S. Basket with vacuum gauge and control valve

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,