

Product Name :
Ideal Delivery Aid for Obstetricians : Safe for Mother / baby**Product Code :**
Gynecological0012**Description :**

Ideal Delivery Aid for Obstetricians : Safe for Mother / baby Housing: ABS Plastic Moulded Capacity: -600 mm Hg \pm 10 , 12 Ltrs/ min Pump Type: Piston Pump Jars: 1 x 1 Ltr Polycarbonate Jar Filter: Bacterial Filter Autoclavable/Reusable Tubing: 8 mm ID x 2 mtr. (PVC) Vacuum Gauge: 5 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 90Watt (110V on request) Noise Level: 55 dB A \pm 3 Dimension & Weight: 30 x 14 x 33 cms, 4.5 Kg Overflow Safety: Mechanical / Electronic Acoustic Alarm Venteuse cup: Silicone - 50, 60 & 70mm : minimizes risk of trauma and cephalohematoma of fetal cranium

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,