

**Product Name :**  
Hospital Platelet Agitator**Product Code :**  
NMBBE10007**Description :**

Hospital Platelet Agitator - Platelet Agitator features the unique design Labtop rack with high ventilation capability to allow oxygen transfer to the platelets. Smooth sliding stainless steel trays which can be removed for cleaning as well as enable easy removal or loading of bags in the tray. Fixed speed agitator with heavy duty bearings and quiet operation. When used with a labtop Platelet Incubator the agitator will automatically stop when the door is open.

No. of Plastic Bags: 24, 36, 48

No. Of tray: 4, 6, 8 Hospital Platelet Agitator, Blood Bank Equipments Exporters, Blood Bank Equipments Manufacturers