

Product Name :
Ceiling Shadowless OT Light with Single Dome Dia 550mm
with 4 Reflectors**Product Code :**
NMOTL10018**Description :**

Ceiling Shadowless OT Light with Single Dome Dia 550mm with 4 Reflectors The light has daylight characters and can be adjusted as per the operation requirements. Lights arms are made of Stainless Steel Fitting. The range of the color temp - 4,200 K \pm 300 as natural color tone Detachable handle for simplistic sterilization Parabolic reflectors & color correcting filters Technical specifications Dome head : 550 mm Halogen Bulbs : 24V x 70W x 4 nos. Light Intensity at 1 mt : 50,000 Lux Height Adjustment : 914 mm Parabolic Reflector : 180 mm Color Temperature : 4200 K \pm 300 Temp rise in field : 2^o - 3^o C Control Unit : Transformer

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [SITE](#)