

**Product Name :**  
Led Shadowless Operation Lamp**Product Code :**  
NMOTL10051**Description :**

Led Shadowless Operation Lamp Spot center illumination : 120000LUX Color temperature : 4300±500K Color indication index (Ra) : 95 Power of LED bead : N64x1W Spot diameter d50 : 130mm-350mm Illumination depth : 120cm Illumination-adjustment range : 10%~100% Ultraviolet irradiance with wavelength under 400nm : 0.002W/m<sup>2</sup> Focusing control : manual Power supply : 220V, 50Hz /110V, 60Hz LED average service life : 60000hrs Lowest height of installation : 2800mm Available in ceiling, floor and wall configurations

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,