

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ceiling Shadowless OT Light Single Dome (5 Bulbs)

Product Code :
NMOTL10055



Description :

Ceiling Shadowless OT Light Single Dome (5 Bulbs) Illumination Intensity : 30000-50000 Lux Color Temperature : 3500-5000 K Spot Diameter : 100-200mm Height Adjustment : 260mm Lamphead Diameter : 500mm Power Supply : 220V / 50Hz Power Input : 150 VA

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: SITE