

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ceiling Mounted Shadowless Operation Lamp With Single Reflector (With Focus Control)

Product Code :
NMOTL10060



Description :

Ceiling Mounted Shadowless Operation Lamp With Single Reflector (With Focus Control) Illumination Intensity : 30000-60000 Lux Color Temperature : 3000-4500 K Spot Diameter : 80-150mm Height Adjustment : 700mm Lamphead Diameter : 195mm Power Supply : 200V / 50Hz Power Input : 75 VA Specification of Halogen Bulb : 50Watt / 24 Volt Brightness Control : Digital Control

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: SITE