



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Ceiling Mounted Led Surgical Light Dual Arm - 2 Lamp

Product Code :
NMOTL10061



Description :

Ceiling Mounted Led Surgical Light Dual Arm - 2 Lamp Diameter of lamp head : Main lamp, 700mm, Subsidiary, 600mm Illumination : Main lamp, 100000 lx ~ 160000 lx, Subsidiary, 80000 lx ~ 160000 lx Color temperature : 3000k ~ 6000k Color rendering index : 92 ~ 100 Diameter of Spot : D10 Main lamp 150mm, Subsidiary 130mm : D50 Main lamp 85mm, Subsidiary 75mm Depth of illumination : = 600mm Adjustment of Height : 800mm Brightness control : 9 steps Total radiation illumination : = 500W/M2 Power supply : A.C100 ~ 220V 50 HZ Power consumption : 120VA Power of bulb : 3 WLED Number of LED Bulbs : Main lamp, 38, Subsidiary 26 Minimum installation height : 3m



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,