



Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Ceiling Shadowless Operation Lamp 2 Arms Double Reflector

**Product Code :**  
NMOTL10064



**Description :**

Ceiling Shadowless Operation Lamp 2 Arms Double Reflector Illumination Intensity : 80000-150000 Lux Color Temperature : 3500-5000 K Spot Diameter : 100-250mm Illumination Depth : 600mm Height Adjustment : 800mm Lamphead Diameter : 700/600mm Brightness Control : 9-level dimmable Temperature Rise of Field : ? 2 (°C) Power Supply : 220V / 50Hz Power Input : 400 VA Specification of Bulb : 150Watt / 24 Volt Minimum Height for Installation : 2750mm Automatic Switchover Between the Key and Backup Bulbs : Yes



**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)