

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Shadowless Operation Lamp With Single Reflector Stand  
Model With Castors (With Focus Control)

**Product Code :**  
NMOTL10071



**Description :**

Shadowless Operation Lamp With Single Reflector Stand Model With Castors (With Focus Control) Illumination Intensity : 30000-60000 Lux Color Temperature : 3000-4500 K Spot Diameter : 80-150mm Height Adjustment : 700mm Lamphead Diameter : 195mm Power Supply : 200V / 50Hz Power Input : 75 VA Specification of Halogen Bulb : 50Watt / 24 Volt Brightness Control : Digital Control

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,