

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
Mobile Shadowless Operation Lamp Single Dome Stand  
Model With Castors & With Battery Backup (4 Bulbs)

**Product Code :**  
NMOTL10074



**Description :**

Mobile Shadowless Operation Lamp Single Dome Stand Model With Castors & With Battery Backup (4 Bulbs)  
Illumination Intensity : 40000 Lux Color Temperature : 3500-5000 K Spot Diameter : 100-200mm Height  
Adjustment : 400mm Lamphead Diameter : 500mm Power Supply : 220V / 50Hz Power Input : 120 VA With  
battery backup

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)