

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Halogen Examination Lamp

Product Code :
NMOTL10076



Description :

Halogen Examination Lamp Illumination Intensity : 20000~30000 Lux Color Temperature : 3000-4000 K Spot Diameter : 80-150 mm Height Adjustment : 600 mm Power Supply : 220V / 50Hz Power Input : 75 VA Specification of Bulb : 50 Watt / 24 Volt

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,