

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
LED Shadowless OT Light Single Reflector Stand Floor
Model With Castors

Product Code :
NMOTL10077



Description :

LED Shadowless OT Light Single Reflector Stand Floor Model With Castors Illumination Intensity : 40000-65000 Lux Color Temperature : 3800-5000 K Color rendering index : 85 (Ra) Spot diameter : 80-160mm Height Adjustment : 700mm Temperature Rise of Field : 0 (°C) Specification of Bulb : LED Bulb 5x3W Power Input : 15 VA Power Supply : 220 V / 50 Hz. AC Brightness Control : Available

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)