

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
LED Shadowless OT Light Single Reflector (Wall Mounted)

**Product Code :**  
NMOTL10078



**Description :**

LED Shadowless OT Light Single Reflector (Wall Mounted) Illumination Intensity : 40000-65000 Lux Color Temperature : 3800-5000 K Color rendering index : 85 (Ra) Spot diameter : 80-160mm Height Adjustment : 700mm Temperature Rise of Field : 0 (°C ) Specification of Bulb : LED Bulb 5x3W Power Input : 15 VA Power Supply : 220 V / 50 Hz. AC Brightness Control : Available

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)