

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
LED Surgical Light Wall Mounted (for Minor Surgical Procedures)

Product Code :
NMOTL10080



Description :

LED Surgical Light Wall Mounted (for Minor Surgical Procedures) Illumination Intensity : 40000 Lux Color Temperature : 6000-6500 K Color rendering index : 75 (Ra) Spot diameter : 150mm Light Field Diameter : D 30 cm = 18 cm Temperature Rise of Field : 0 (°C) Power Supply : 220 V / 50 Hz. AC Extention : 75 cm goose neck Specification of Bulb : LED 6 bulbs (high power) Brightness Control : Available

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)