

**NAUGRA**<sup>®</sup>

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
LED Examination Light (Wall Mounted)

**Product Code :**  
NMOTL10082



**Description :**

LED Examination Light (Wall Mounted) Illumination Intensity : 12000 Lux Color Temperature : 6000-6500 K Spot Diameter : 80 - 250mm Power Supply : 220 V / 50 Hz. AC Extension : 75 cm goose neck Specification of Bulb : LED 1x3W Brightness Control : Available

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)