

**Product Name :**  
LED Examination Light With Stand Floor Model & With  
Microcomputer Controlled 8 Step Brightness Adjustment

**Product Code :**  
NMOTL10083



**Description :**

LED Examination Light With Stand Floor Model & With Microcomputer Controlled 8 Step Brightness Adjustment  
Standard Features : Microcomputer control 8-steps brightness adjustment Technical specifications Color  
Temperature : 3000-5500 K Illumination : ? 20000 Lux Temperature Rise of Field : ?8(°C) Power Supply : 220 V  
/ 50 Hz. AC Power (W) : ?12 Type of bulb : Super bright LED Luminotron 3V/0.5Wx19 Spot Diameter : 80-120  
mm Range : 300-500mm from lamp adapter Brightness Control : 8 step digital control

**NAUGRA**<sup>®</sup>

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)