

Product Name :
LED Examination Light With Stand Floor Model

Product Code :
NMOTL10085



Description :

LED Examination Light With Stand Floor Model Color Temperature : 3000-5500 K Illumination : ? 12000 Lux
Temperature Rise of Field : ? 8 (°C) Power Supply : 220 V / 50 Hz. AC Power : ? 3.5 (W) Type of bulb : Super
bright LED Luminotron 3V/0.5Wx8 Spot Diameter : 80-120mm Range : 300-500mm from lamp adapter
Brightness Control : Switch ON/OFF

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,