

**Product Name :**  
LED Examination Light With Stand Floor Model

**Product Code :**  
NMOTL10085



**Description :**

LED Examination Light With Stand Floor Model Color Temperature : 3000-5500 K Illumination : ? 12000 Lux  
Temperature Rise of Field : ? 8 (°C) Power Supply : 220 V / 50 Hz. AC Power : ? 3.5 (W) Type of bulb : Super  
bright LED Luminotron 3V/0.5Wx8 Spot Diameter : 80-120mm Range : 300-500mm from lamp adapter  
Brightness Control : Switch ON/OFF

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)