

NAUGRA[®]

Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
LED Examination Light Table Top Model

Product Code :
NMOTL10087



Description :

LED Examination Light Table Top Model Technical specifications Color Temperature : 3000-5500 K Illumination : ? 12000 Lux Temperature Rise of Field : ? 8 (°C) Power Supply : A.C 220V 50Hz Power : ? 3.5 (W) Type of bulb : Super bright LED Luminotron 3V/0.5Wx8 Spot Diameter : 80-120mm Range : 300-500mm from lamp adapter Brightness Control : Digital

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website: [SITE](#)