

**NAUGRA<sup>®</sup>**

Email : sales@naugramedical.com  
Phone: 9896600003

**Product Name :**  
LED Examination Light

**Product Code :**  
NMOTL10089



**Description :**

LED Examination Light Technical specifications Color Temperature : 3000-5500 K Illumination : ? 12000 Lux  
Temperature Rise of Field : ? 8 (°C) Power Supply : A.C 220V 50Hz Power : ? 3.5 (W) Type of bulb : Super  
bright LED Luminotron 3V/0.5Wx8 Spot Diameter : 80-120mm Range : 300-500mm from lamp adapter  
Brightness Control : Digital

**NAUGRA<sup>®</sup>**

**Medical Hospital Equipments Manufacturers**

Website: [SITE](#)