

Product Name :
Suction Unit - INSTAD AC/DC**Product Code :**
NMSUC20HV4**Description :**

Suction Unit - INSTAD AC/DC Housing: ABS Plastic Moulded & 3 Stage LED Indicator for Battery Status
Capacity: -600 mm Hg \pm 10 , 21 Ltrs/ min Pump Type: Piston Pump Jars: 1 x 1 Ltr Polycarbonate Jar Filter:
Bacterial Filter Autoclavable/Reusable Tubing: 8 mm ID x 2 mtr. (Silicon) Vacuum Gauge: 5 cm dia 0-760 mm
Hg calibration. Power: 220 V AC, 50 / 60 Hz, 70Watt , 12 V DC (Internal Rechargeable Battery 12V 4.5Amp)
Battery: Lead Acid (Ni-cd Optional) Noise Level: 55 dB A \pm 3 Dimension & Weight: 39 x 22 x 22 cms, 6 Kg.
Overflow Safety: Mechanical Accessories: Car Charger & Power Adaptor

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,