

Product Name :
Suction Unit – Instad DC**Product Code :**
NMSUC20HV5**Description :**

Suction Unit – Instad DC Housing: ABS Plastic Moulded Capacity: -600 mm Hg \pm 10 , 21 Ltrs/ min Pump Type: Piston Pump Jars: 1 x 1 Ltr Polycarbonate Jar Filter: Bacterial Filter Autoclavable/Reusable Tubing: 8 mm ID x 2 mtr. (Silicon) Vacuum Gauge: 5 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 70Watt , 12 V DC (No Internal Battery) Noise Level: 55 dB A \pm 3 Dimension & Weight: 39 x 22 x 22 cms, 4 Kg. Overflow Safety: Mechanical Accessories: Car Charger & Power Adaptor

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE