

Product Name :
Suction Unit - HIVAC**Product Code :**
NMSUC10HV4**Description :**

Suction Unit - HIVAC Housing: MS Powder Coated on four castors Capacity: -700 mm Hg \pm 10, 50 Ltrs/ min
Pump Type: Double Rotary Vane Pump Jars: 2 x 2 Ltrs Polycarbonate Jar Filters: Trap Filter Tubing: 10 mm ID x
2 mtr. (PVC) Vacuum Gauge: 10 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 380Watt (110V
on request) Noise Level: 50 dB A \pm 3 Dimension & Weight: 49 x 36 x 70 cms, 37 Kg. Overflow Safety:
Mechanical

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,