

**Product Name :**  
Suction Unit - HIVAC**Product Code :**  
NMSUC10HV4**Description :**

Suction Unit - HIVAC Housing: MS Powder Coated on four castors Capacity: -700 mm Hg  $\pm$  10, 50 Ltrs/ min  
Pump Type: Double Rotary Vane Pump Jars: 2 x 2 Ltrs Polycarbonate Jar Filters: Trap Filter Tubing: 10 mm ID x  
2 mtr. (PVC) Vacuum Gauge: 10 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 380Watt (110V  
on request) Noise Level: 50 dB A  $\pm$  3 Dimension & Weight: 49 x 36 x 70 cms, 37 Kg. Overflow Safety:  
Mechanical

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)